

产品名称: Alpha 1 adaptin Rabbit Monoclonal Antibody
产品货号: AMRe84435



产品概述 (Summary)

产品名称 (Production Name)	Alpha 1 adaptin Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF,ICC,FC
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide,0.05% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	Alpha 1 adaptin
别名 (Alternative Names)	ADTAA; CLAPA1;;AP 2 complex alpha 1
基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	O95782.A synthesized peptide derived from human AP 2 complex alpha 1

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:2000,IHC 1:100-1:200,ICC/IF 1:50-1:200,ICC 1:50-1:200,FC 1:20-1:100
蛋白分子量 (Molecular Weight)	108 kDa

研究背景 (Background)


产品名称: Alpha 1 adaptin Rabbit Monoclonal Antibody
产品货号: AMRe84435




Component of the adaptor protein complex 2 (AP-2). Adaptor protein complexes function in protein transport via transport vesicles in different membrane traffic pathways. Adaptor protein complexes are vesicle coat components and appear to be involved in cargo selection and vesicle formation.

研究领域 (Research Area)

图片 (Image Data)



All lanes use the Antibody at 1:2K dilution for 1 hour at room temperature.



Western blot analysis of Alpha1-adaptin expression in A431 cell lysate.

注意事项 (Note)

For research use only .