

产品名称: NUP153 Rabbit Monoclonal Antibody
产品货号: AMRe84437



产品概述 (Summary)

产品名称 (Production Name)	NUP153 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB, ICC, FC
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide, 0.05% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	NUP153
别名 (Alternative Names)	HNUP153; Nup153;;NUP153
基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	P49790.A synthesized peptide derived from human NUP153

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:2000, ICC 1:50-1:200, FC 1:20-1:100
蛋白分子量 (Molecular Weight)	154 kDa

研究背景 (Background)

Possible DNA-binding subunit of the nuclear pore complex (NPC). The repeat-containing domain may be involved in


产品名称: NUP153 Rabbit Monoclonal Antibody
产品货号: AMRe84437



anchoring components of the pore complex to the pore membrane.

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of Nup153 expression in K562 cell lysate.

注意事项 (Note)

For research use only .