

产品名称: Phospho-PKC(S729) Rabbit Monoclonal Antibody
产品货号: AMRe84443



产品概述 (Summary)

产品名称 (Production Name)	Phospho-PKC(S729) Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB
种属反应性 (Reactivity)	Human,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Phosphorylated
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide, 0.05% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	Phospho-PKC(S729)
别名 (Alternative Names)	PKC alpha; PKC beta; PKC delta; PKC epsilon; PKC gamma; PKC zeta; PKC2; PKCA; PKCB; PKCD; PKCE; PKCG; PRKCA; PRKCB; PRKCD; PRKCE; PRKCG; PRKCZ;;p-PKC epsilon (S729)
基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	Q02156(PKCepsilon),P05771(PKCbeta),P17252(PKCalpha).A synthesized peptide derived from human PKC epsilon around the phosphorylation site of S729

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:2000
----------------------	------------------

产品名称: Phospho-PKC(S729) Rabbit Monoclonal Antibody
产品货号: AMRe84443




蛋白分子量 (Molecular Weight) 84 kDa

研究背景 (Background)

Calcium-independent, phospholipid- and diacylglycerol (DAG)-dependent serine/threonine-protein kinase that plays essential roles in the regulation of multiple cellular processes linked to cytoskeletal proteins, such as cell adhesion, motility, migration and cell cycle, functions in neuron growth and ion channel regulation, and is involved in immune response, cancer cell invasion and regulation of apoptosis.

研究领域 (Research Area)

图片 (Image Data)



All lanes use the Antibody at 1:2K dilution for 1 hour at room temperature.

注意事项 (Note)

For research use only .