

产品名称: ESRRG Rabbit Monoclonal Antibody
产品货号: AMRe84450



产品概述 (Summary)

产品名称 (Production Name)	ESRRG Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide,0.05% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	ESRRG
别名 (Alternative Names)	ERR 3; ERR gamma 2; ERR3; ERRG2; ERRgamma; ESRRG; NR3B3;;Estrogen related receptor gamma
基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	P62508.A synthesized peptide derived from human Estrogen related receptor gamma

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:2000
蛋白分子量 (Molecular Weight)	51 kDa

产品名称: ESRRG Rabbit Monoclonal Antibody
产品货号: AMRe84450




研究背景 (Background)


Orphan receptor that acts as a transcription activator in the absence of bound ligand. Binds specifically to an estrogen response element and activates reporter genes controlled by estrogen response elements.

研究领域 (Research Area)

图片 (Image Data)



All lanes use the Antibody at 1:1K dilution for 1 hour at room temperature.



Western blot analysis of Estrogen Related Receptor gamma expression in Human fetal heart lysate.

注意事项 (Note)

For research use only .