

产品名称: Laminin gamma 1 Rabbit Monoclonal Antibody  
产品货号: AMRe84452



## 产品概述 (Summary)

|                        |                                            |
|------------------------|--------------------------------------------|
| 产品名称 (Production Name) | Laminin gamma 1 Rabbit Monoclonal Antibody |
| 描述 (Description)       | Recombinant rabbit monoclonal antibody     |
| 宿主 (Host)              | Rabbit                                     |
| 应用 (Application)       | WB,FC                                      |
| 种属反应性 (Reactivity)     | Human                                      |

## 产品性能 (Performance)

|                     |                                                                                              |
|---------------------|----------------------------------------------------------------------------------------------|
| 偶联物 (Conjugation)   | Unconjugated                                                                                 |
| 修饰 (Modification)   | Unmodified                                                                                   |
| 同种型 (Isotype)       | IgG                                                                                          |
| 克隆 (Clonality)      | Monoclonal                                                                                   |
| 形式 (Form)           | Liquid                                                                                       |
| 存放说明 (Storage)      | Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.     |
| 储存溶液 (Buffer)       | Purified antibody in PBS with 0.05% sodium azide, 0.05% protective protein and 50% glycerol. |
| 纯化方式 (Purification) | Affinity Purification                                                                        |

## 免疫原信息 (Immunogen)

|                        |                                                                                                                                                                                                |
|------------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| 基因名 (Gene Name)        | Laminin gamma 1<br>LAMB2; lamc1; Laminin B2 chain; Laminin subunit gamma-1; Laminin-1 subunit gamma; Laminin-10 subunit gamma; Laminin-11 subunit gamma;                                       |
| 别名 (Alternative Names) | Laminin-2 subunit gamma; Laminin-3 subunit gamma; Laminin-4 subunit gamma; Laminin-6 subunit gamma; Laminin-7 subunit gamma; Laminin-8 subunit gamma; Laminin-9 subunit gamma;;Laminin gamma 1 |
| 基因 ID (Gene ID)        |                                                                                                                                                                                                |
| 蛋白 ID (SwissProt ID)   | P11047/Q5VYE7/Q6NVY8.A synthesized peptide derived from human Laminin gamma 1                                                                                                                  |

## 产品应用 (Application)

产品名称: Laminin gamma 1 Rabbit Monoclonal Antibody  
产品货号: AMRe84452



稀释比 (Dilution Ratio)

WB 1:1000-1:2000, FC 1:20-1:100

蛋白分子量 (Molecular Weight)


Calculated MW: 178 kDa ; Observed MW: 200-280 kDa

## 研究背景 (Background)


Binding to cells via a high affinity receptor, laminin is thought to mediate the attachment, migration and organization of cells into tissues during embryonic development by interacting with other extracellular matrix components.

## 研究领域 (Research Area)

## 图片 (Image Data)



Western blot analysis of Laminin gamma 1 expression in A431 cell lysate.



All lanes use the Antibody at 1:1k dilution for 1 hour at room temperature.

## 注意事项 (Note)

For research use only .