

产品名称: Cytochrome P450 4A Rabbit Monoclonal Antibody
产品货号: AMRe84455



产品概述 (Summary)

产品名称 (Production Name)	Cytochrome P450 4A Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,IP
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide,0.05% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	Cytochrome P450 4A
别名 (Alternative Names)	CYP4A; Cyp4a1; Cyp4a10; CYP4A11; Cyp4a14; Cyp4a3; CYP4A7; CYP4A11; CYP4A11;;Cytochrome P450 4A11
基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	Q02928.A synthesized peptide derived from human Cytochrome P450 4A11

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:2000,IHC 1:100-1:200,IP 1:20-1:50
蛋白分子量 (Molecular Weight)	Calculated MW: 59 kDa ; Observed MW: 50 kDa

研究背景 (Background)


产品名称: Cytochrome P450 4A Rabbit Monoclonal Antibody
产品货号: AMRe84455



Catalyzes the omega- and (omega-1)-hydroxylation of various fatty acids such as laurate, myristate and palmitate. Has little activity toward prostaglandins A1 and E1. Oxidizes arachidonic acid to 20-hydroxyecosatetraenoic acid (20-HETE).

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of Cytochrome P450 4A expression in Human fetal kidney lysate.

注意事项 (Note)

For research use only .