

产品名称: COL11A1 Rabbit Monoclonal Antibody
产品货号: AMRe84457



产品概述 (Summary)

产品名称 (Production Name)	COL11A1 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide, 0.05% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	COL11A1
别名 (Alternative Names)	COL11A1; COLL6; Collagen alpha 1; collagen XI alpha 1; STL2; STL3;; Collagen XI alpha 1
基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	P12107.A synthesized peptide derived from human Collagen XI alpha 1

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:2000
蛋白分子量 (Molecular Weight)	Calculated MW: 181 kDa ; Observed MW: 150 kDa

研究背景 (Background)


产品名称: COL11A1 Rabbit Monoclonal Antibody
产品货号: AMRe84457



May play an important role in fibrillogenesis by controlling lateral growth of collagen II fibrils.

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of COL11A1 expression in K562 cell lysate.

注意事项 (Note)

For research use only .