

产品名称: FGF19 Rabbit Monoclonal Antibody
产品货号: AMRe84461



产品概述 (Summary)

产品名称 (Production Name)	FGF19 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide, 0.05% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	FGF19
别名 (Alternative Names)	FGF15; FGF19; Fibroblast growth factor 15; Fibroblast growth factor 19;;FGF19
基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	O95750.A synthesized peptide derived from human FGF19

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:2000
蛋白分子量 (Molecular Weight)	Calculated MW: 24 kDa ; Observed MW: 36 kDa

研究背景 (Background)

Involved in the suppression of bile acid biosynthesis through down-regulation of CYP7A1 expression, following positive


产品名称: FGF19 Rabbit Monoclonal Antibody
产品货号: AMRe84461



regulation of the JNK and ERK1/2 cascades. Stimulates glucose uptake in adipocytes. Activity requires the presence of KLB and FGFR4.

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of FGF19 expression in Human fetal liver lysate.

注意事项 (Note)

For research use only .