

产品名称: CNTP2 Rabbit Monoclonal Antibody
产品货号: AMRe84471



产品概述 (Summary)

产品名称 (Production Name)	CNTP2 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide,0.05% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	CNTP2
别名 (Alternative Names)	AUTS15; CDFE; CNTNAP2; CNTP2; NRXN4; PTHSL1;;CNTNAP2
基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	Q9UHC6.A synthesized peptide derived from human CNTNAP2

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:2000,IHC 1:100-1:200
蛋白分子量 (Molecular Weight)	Calculated MW: 148 kDa ; Observed MW: 160 kDa

研究背景 (Background)

Required for gap junction formation (Probable). Required, with CNTNAP1, for radial and longitudinal organization of


产品名称: CNTP2 Rabbit Monoclonal Antibody
产品货号: AMRe84471




myelinated axons. Plays a role in the formation of functional distinct domains critical for saltatory conduction of nerve impulses in myelinated nerve fibers. Demarcates the juxtaparanodal region of the axo-glial junction.

研究领域 (Research Area)


图片 (Image Data)



All lanes use the Antibody at 1:2K dilution for 1 hour at room temperature.



All lanes use the Antibody at 1:1K dilution for 1 hour at room temperature.



All lanes use the Antibody at 1:1K dilution for 1 hour at room temperature.

产品名称: CNTP2 Rabbit Monoclonal Antibody
产品货号: AMRe84471



注意事项 (Note)

For research use only .