

产品名称: BAF57/SMARCE1 Rabbit Monoclonal Antibody  
产品货号: AMRe84474



## 产品概述 (Summary)

|                        |                                          |
|------------------------|------------------------------------------|
| 产品名称 (Production Name) | BAF57/SMARCE1 Rabbit Monoclonal Antibody |
| 描述 (Description)       | Recombinant rabbit monoclonal antibody   |
| 宿主 (Host)              | Rabbit                                   |
| 应用 (Application)       | WB,IHC,ICC/IF,ICC,IP                     |
| 种属反应性 (Reactivity)     | Human,Mouse,Rat                          |

## 产品性能 (Performance)

|                     |                                                                                             |
|---------------------|---------------------------------------------------------------------------------------------|
| 偶联物 (Conjugation)   | Unconjugated                                                                                |
| 修饰 (Modification)   | Unmodified                                                                                  |
| 同种型 (Isotype)       | IgG                                                                                         |
| 克隆 (Clonality)      | Monoclonal                                                                                  |
| 形式 (Form)           | Liquid                                                                                      |
| 存放说明 (Storage)      | Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.    |
| 储存溶液 (Buffer)       | Purified antibody in PBS with 0.05% sodium azide,0.05% protective protein and 50% glycerol. |
| 纯化方式 (Purification) | Affinity Purification                                                                       |

## 免疫原信息 (Immunogen)

|                        |                                                         |
|------------------------|---------------------------------------------------------|
| 基因名 (Gene Name)        | BAF57/SMARCE1                                           |
| 别名 (Alternative Names) | BAF57; SMARCE1;;SMARCE1                                 |
| 基因 ID (Gene ID)        |                                                         |
| 蛋白 ID (SwissProt ID)   | Q969G3.A synthesized peptide derived from human SMARCE1 |

## 产品应用 (Application)

|                          |                                                                                |
|--------------------------|--------------------------------------------------------------------------------|
| 稀释比 (Dilution Ratio)     | WB 1:1000-1:2000,IHC 1:100-1:200,ICC/IF 1:50-1:200,ICC 1:50-1:200,IP 1:20-1:50 |
| 蛋白分子量 (Molecular Weight) | Calculated MW: 47 kDa ; Observed MW: 57 kDa                                    |

## 研究背景 (Background)


**产品名称: BAF57/SMARCE1 Rabbit Monoclonal Antibody**  
**产品货号: AMRe84474**

---

Involved in transcriptional activation and repression of select genes by chromatin remodeling (alteration of DNA-nucleosome topology). Component of SWI/SNF chromatin remodeling complexes that carry out key enzymatic activities, changing chromatin structure by altering DNA-histone contacts within a nucleosome in an ATP-dependent manner. Belongs to the neural progenitors-specific chromatin remodeling complex (npBAF complex) and the neuron-specific chromatin remodeling complex (nBAF complex).

## 研究领域 (Research Area)

## 图片 (Image Data)



Western blot analysis of BAF57/SMARCE1 expression in Raji cell lysate.

## 注意事项 (Note)

For research use only .