

产品名称: ADK Rabbit Monoclonal Antibody
产品货号: AMRe84475



产品概述 (Summary)

产品名称 (Production Name)	ADK Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB
种属反应性 (Reactivity)	Human,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide, 0.05% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	ADK
别名 (Alternative Names)	Adenosine kinase; adk; AK;;ADK
基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	P55263.A synthesized peptide derived from human ADK

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:2000
蛋白分子量 (Molecular Weight)	Calculated MW: 41 kDa ; Observed MW: 46 kDa

研究背景 (Background)

ATP dependent phosphorylation of adenosine and other related nucleoside analogs to monophosphate derivatives. Serves


产品名称: ADK Rabbit Monoclonal Antibody
产品货号: AMRe84475




as a potential regulator of concentrations of extracellular adenosine and intracellular adenine nucleotides.

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of ADK expression in HepG2 cell lysate.



All lanes use the Antibody at 1:1K dilution for 1 hour at room temperature.

注意事项 (Note)

For research use only .