

产品名称: HCadherin Rabbit Monoclonal Antibody
产品货号: AMRe84490



产品概述 (Summary)

产品名称 (Production Name)	HCadherin Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide,0.05% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	HCadherin
别名 (Alternative Names)	CDH13; CDHH; H-cadherin; Heart cadherin; P105; T cad; T Cadherin; Truncated cadherin;;Cadherin 13
基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	P55290.A synthesized peptide derived from human Cadherin 13

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:2000,IHC 1:100-1:200
蛋白分子量 (Molecular Weight)	Calculated MW: 78 kDa ; Observed MW: 105,130 kDa

研究背景 (Background)


产品名称: HCadherin Rabbit Monoclonal Antibody
产品货号: AMRe84490




Cadherins are calcium dependent cell adhesion proteins. They preferentially interact with themselves in a homophilic manner in connecting cells; cadherins may thus contribute to the sorting of heterogeneous cell types. May act as a negative regulator of neural cell growth.

研究领域 (Research Area)


图片 (Image Data)



Western blot analysis of H Cadherin expression in Human fetal heart lysate.




Immunohistochemical analysis of paraffin-embedded Mouse skeletal muscle, using the Antibody at 1:100 dilution.



Immunohistochemical analysis of paraffin-embedded Human cervical cancer, using the Antibody at 1:100 dilution.

产品名称: HCadherin Rabbit Monoclonal Antibody
产品货号: AMRe84490



Immunohistochemical analysis of paraffin-embedded Human squamous cell carcinoma , using the Antibody at 1:100 dilution.

注意事项 (Note)

For research use only .