

产品名称: PDE2A Rabbit Monoclonal Antibody
产品货号: AMRe84493



产品概述 (Summary)

产品名称 (Production Name)	PDE2A Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide,0.05% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	PDE2A
别名 (Alternative Names)	CGS PDE; cGSPDE; PDE2A; PDE2A1; PDE2A4; phosphodiesterase 2A;;PDE2A
基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	O00408.A synthesized peptide derived from human PDE2A

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:2000
蛋白分子量 (Molecular Weight)	106 kDa

研究背景 (Background)

cGMP-activated cyclic nucleotide phosphodiesterase with a dual-specificity for the second messengers cAMP and cGMP,


产品名称: PDE2A Rabbit Monoclonal Antibody
产品货号: AMRe84493



which are key regulators of many important physiological processes.

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of PDE2A expression in Human fetal brain lysate.

注意事项 (Note)

For research use only .