

产品名称: Mad2L1 Rabbit Monoclonal Antibody
产品货号: AMRe84499



产品概述 (Summary)

产品名称 (Production Name)	Mad2L1 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide, 0.05% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	Mad2L1
别名 (Alternative Names)	HsMAD2; MAD 2; MAD 2A; MAD2 like 1; Mad2l1;;HsMAD2
基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	Q13257.A synthesized peptide derived from human HsMAD2

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:2000
蛋白分子量 (Molecular Weight)	24 kDa

研究背景 (Background)

Component of the spindle-assembly checkpoint that prevents the onset of anaphase until all chromosomes are properly


产品名称: Mad2L1 Rabbit Monoclonal Antibody
产品货号: AMRe84499



aligned at the metaphase plate. Required for the execution of the mitotic checkpoint which monitors the process of kinetochore-spindle attachment and inhibits the activity of the anaphase promoting complex by sequestering CDC20 until all chromosomes are aligned at the metaphase plate.

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of Mad2L1 expression in A431 cell lysate.

注意事项 (Note)

For research use only .