

产品名称: AP1G1 Rabbit Monoclonal Antibody
产品货号: AMRe84501



产品概述 (Summary)

产品名称 (Production Name)	AP1G1 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide,0.05% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	AP1G1
别名 (Alternative Names)	ADTG; AP1G1; CLAPG1;;gamma 1 Adaptin
基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	O43747.A synthesized peptide derived from human gamma 1 Adaptin

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:2000,IHC 1:100-1:200
蛋白分子量 (Molecular Weight)	Calculated MW: 91 kDa ; Observed MW: 100 kDa

研究背景 (Background)

Subunit of clathrin-associated adaptor protein complex 1 that plays a role in protein sorting in the late-Golgi/trans-Golgi


产品名称: AP1G1 Rabbit Monoclonal Antibody
产品货号: AMRe84501




network (TGN) and/or endosomes. The AP complexes mediate both the recruitment of clathrin to membranes and the recognition of sorting signals within the cytosolic tails of transmembrane cargo molecules.

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of AP1G1 expression in 293T cell lysate.



All lanes use the Antibody at 1:1K dilution for 1 hour at room temperature.

注意事项 (Note)

For research use only .