

产品名称: MBD3 Rabbit Monoclonal Antibody
产品货号: AMRe84502



产品概述 (Summary)

产品名称 (Production Name)	MBD3 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB, ICC, FC, IP
种属反应性 (Reactivity)	Human, Mouse, Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide, 0.05% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	MBD3
别名 (Alternative Names)	Mbd3;;MBD3
基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	O95983.A synthesized peptide derived from human MBD3

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:2000, ICC 1:50-1:200, FC 1:20-1:100, IP 1:20-1:50
蛋白分子量 (Molecular Weight)	33 kDa

研究背景 (Background)

Acts as a component of the histone deacetylase NuRD complex which participates in the remodeling of chromatin. Acts as


产品名称: MBD3 Rabbit Monoclonal Antibody
产品货号: AMRe84502




transcriptional repressor and plays a role in gene silencing.

研究领域 (Research Area)


图片 (Image Data)



Western blot analysis of MBD3 expression in 293T cell lysate.



All lanes use the Antibody at 1:2K dilution for 1 hour at room temperature.



All lanes use the Antibody at 1:2K dilution for 1 hour at room temperature.

产品名称: MBD3 Rabbit Monoclonal Antibody
产品货号: AMRe84502



注意事项 (Note)

For research use only .