

产品名称: Myogenin Rabbit Monoclonal Antibody
产品货号: AMRe84507



产品概述 (Summary)

产品名称 (Production Name)	Myogenin Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,FC
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide,0.05% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	Myogenin
别名 (Alternative Names)	bHLHc3; cb553; MYF4; MYOG; Myogenic factor 4;;Myogenin
基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	P15173.A synthesized peptide derived from human Myogenin

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:2000,IHC 1:100-1:200,FC 1:20-1:100
蛋白分子量 (Molecular Weight)	Calculated MW: 25 kDa ; Observed MW: 36 kDa

研究背景 (Background)

Acts as a transcriptional activator that promotes transcription of muscle-specific target genes and plays a role in muscle


产品名称: Myogenin Rabbit Monoclonal Antibody
产品货号: AMRe84507



differentiation, cell cycle exit and muscle atrophy. Essential for the development of functional embryonic skeletal fiber muscle differentiation.

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of Myogenin expression in C2C12 cell lysate.

注意事项 (Note)

For research use only .