

产品名称: SCGN Rabbit Monoclonal Antibody
产品货号: AMRe84509



产品概述 (Summary)

产品名称 (Production Name)	SCGN Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,IP
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide,0.05% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	SCGN
别名 (Alternative Names)	CALBL; Scgn; SECRET; Secretagogin; SEGN; Setagin;;Secretagogin
基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	O76038.A synthesized peptide derived from human Secretagogin

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:2000,IHC 1:100-1:200,IP 1:20-1:50
蛋白分子量 (Molecular Weight)	32 kDa

研究背景 (Background)

The encoded protein is a secreted calcium-binding protein which is found in the cytoplasm. It is related to calbindin D-28K


产品名称: SCGN Rabbit Monoclonal Antibody
产品货号: AMRe84509




and calretinin. This protein is thought to be involved in KCL-stimulated calcium flux and cell proliferation.

研究领域 (Research Area)


图片 (Image Data)



All lanes use the Antibody at 1:2K dilution for 1 hour at room temperature.



All lanes use the Antibody at 1:2K dilution for 1 hour at room temperature.



All lanes use the Antibody at 1:2K dilution for 1 hour at room temperature.

产品名称: SCGN Rabbit Monoclonal Antibody
产品货号: AMRe84509



注意事项 (Note)

For research use only .