

产品名称: TBX21 Rabbit Monoclonal Antibody
产品货号: AMRe84511



产品概述 (Summary)

产品名称 (Production Name)	TBX21 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	IHC,FC
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide, 0.05% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	TBX21
别名 (Alternative Names)	T bet; T box 21; T PET; TBET; TBLYM; Tbx21; TPET; Transcription factor TBLYM;;TBX21
基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	Q9UL17.A synthesized peptide derived from human TBX21

产品应用 (Application)

稀释比 (Dilution Ratio)	IHC 1:100-1:200, FC 1:20-1:100
蛋白分子量 (Molecular Weight)	58 kDa

研究背景 (Background)


产品名称: TBX21 Rabbit Monoclonal Antibody
产品货号: AMRe84511



Lineage-defining transcription factor which initiates Th1 lineage development from naive Th precursor cells both by activating Th1 genetic programs and by repressing the opposing Th2 and Th17 genetic programs. Activates transcription of a set of genes important for Th1 cell function, including those encoding IFN-gamma and the chemokine receptor CXCR3. Induces permissive chromatin accessibility and CpG methylation in IFNG.

研究领域 (Research Area)

图片 (Image Data)



Immunohistochemical analysis of paraffin-embedded human spleen, using TBX21 Antibody.

注意事项 (Note)

For research use only .