

产品名称: TFE3 Rabbit Monoclonal Antibody
产品货号: AMRe84512



产品概述 (Summary)

产品名称 (Production Name)	TFE3 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF,ICC,FC
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide,0.05% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	TFE3
别名 (Alternative Names)	bHLHe33; RCCP2; RCCX1; Tcf3; Tfe3; TFEA;;TFE3
基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	P19532.A synthesized peptide derived from human TFE3

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:2000,IHC 1:100-1:200,ICC/IF 1:50-1:200,ICC 1:50-1:200,FC 1:20-1:100
蛋白分子量 (Molecular Weight)	62 kDa

研究背景 (Background)


产品名称: TFE3 Rabbit Monoclonal Antibody
产品货号: AMRe84512




Transcription factor that acts as a master regulator of lysosomal biogenesis and immune response . Specifically recognizes and binds E-box sequences (5'-CANNTG-3'); efficient DNA-binding requires dimerization with itself or with another Mit/TFE family member such as TFEB or MITF.

研究领域 (Research Area)


图片 (Image Data)



All lanes use the Antibody at 1:1K dilution for 1 hour at room temperature.




Immunofluorescent analysis using the Antibody at 1:150 dilution.



Immunofluorescent analysis using the Antibody at 1:50 dilution.

产品名称: TFE3 Rabbit Monoclonal Antibody
产品货号: AMRe84512



Immunohistochemical analysis of paraffin-embedded human bladder, using TFE3 Antibody.

注意事项 (Note)

For research use only .