

产品名称: BCL2L15 Rabbit Monoclonal Antibody
产品货号: AMRe84520



产品概述 (Summary)

产品名称 (Production Name)	BCL2L15 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF,ICC,FC
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide,0.05% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	BCL2L15
别名 (Alternative Names)	Bcl-2-like protein 15; Bcl2-L-15; BCL2-like 15; Bcl2l15; Bfk; Gm566; Pro apoptotic Bcl 2 protein; ;Bcl 2 L15
基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	Q5TBC7.A synthesized peptide derived from human Bcl 2 L15

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:2000,IHC 1:100-1:200,ICC/IF 1:50-1:200,ICC 1:50-1:200,FC 1:20-1:100
蛋白分子量 (Molecular Weight)	Calculated MW: 18 kDa ; Observed MW: 52 kDa

产品名称: BCL2L15 Rabbit Monoclonal Antibody
产品货号: AMRe84520




研究背景 (Background)


Predicted to be involved in apoptotic process and regulation of apoptotic process. Predicted to be active in cytosol and nucleus.

研究领域 (Research Area)


图片 (Image Data)



All lanes use the Antibody at 1:1K dilution for 1 hour at room temperature.



All lanes use the Antibody at 1:1K dilution for 1 hour at room temperature.



Western blot analysis of BCL2L15 in His-tagged BCL2L15 cell lysate.

产品名称: BCL2L15 Rabbit Monoclonal Antibody
产品货号: AMRe84520



注意事项 (Note)

For research use only .