

产品名称: IL15RA Rabbit Monoclonal Antibody
产品货号: AMRe84521



产品概述 (Summary)

产品名称 (Production Name)	IL15RA Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB, ICC
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide, 0.05% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	IL15RA CD215; IL 15R alpha; IL-15 receptor subunit alpha; IL-15R-alpha; IL-15RA;
别名 (Alternative Names)	IL15ra; Interleukin 15 receptor alpha; sIL-15 receptor subunit alpha; sIL-15R-alpha; sIL-15RA ; Soluble interleukin 15 receptor subunit alpha;; IL 15 receptor alpha
基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	Q13261.A synthesized peptide derived from human IL 15 receptor alpha

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:2000, ICC 1:50-1:200
蛋白分子量 (Molecular Weight)	Calculated MW: 28 kDa ; Observed MW: 38 kDa

产品名称: IL15RA Rabbit Monoclonal Antibody
产品货号: AMRe84521




研究背景 (Background)

Receptor for interleukin-15. Expression of different isoforms may alter or interfere with signal transduction. Isoform 5, isoform 6, isoform 7 and isoform 8 do not bind IL15. Signal transduction involves STAT3, STAT5, STAT6, JAK2 and SYK.

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of IL15RA expression in Raji cell lysate.

注意事项 (Note)

For research use only .