

产品名称: Phospho-RIP3(S232) Rabbit Monoclonal Antibody
产品货号: AMRe84524



产品概述 (Summary)

产品名称 (Production Name)	Phospho-RIP3(S232) Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB
种属反应性 (Reactivity)	Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Phosphorylated
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide, 0.05% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	Phospho-RIP3(S232) Receptor interacting protein 3; Receptor interacting serine threonine kinase 3;
别名 (Alternative Names)	Receptor interacting serine/threonine protein kinase 3; Receptor-interacting serine/threonine-protein kinase 3; RIP like protein kinase 3; RIP-3; RIP-like protein kinase 3; RIPK3; ;p-Rip3 (S232)
基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	Q9QZL0(mouse)(human:Q9Y572).A synthesized peptide derived from mouse Rip3 around the phosphorylation site of S232

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:2000
----------------------	------------------

产品名称: Phospho-RIP3(S232) Rabbit Monoclonal Antibody
产品货号: AMRe84524




蛋白分子量 (Molecular Weight) Calculated MW: 53 kDa ; Observed MW: 50 kDa

研究背景 (Background)

Serine/threonine-protein kinase that activates necroptosis and apoptosis, two parallel forms of cell death. Necroptosis, a programmed cell death process in response to death-inducing TNF- α family members, is triggered by RIPK3 following activation by ZBP1.

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of Phospho-RIP3 (S232) expression in L929 cell lysate.

注意事项 (Note)

For research use only .