

产品名称: SYVN1/HRD1 Rabbit Monoclonal Antibody
产品货号: AMRe84525



产品概述 (Summary)

产品名称 (Production Name)	SYVN1/HRD1 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB, ICC, FC, IP
种属反应性 (Reactivity)	Human, Mouse, Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide, 0.05% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	SYVN1/HRD1
别名 (Alternative Names)	DER3; E3 ubiquitin-protein ligase synoviolin; HRD1; Synovial apoptosis inhibitor 1; Synoviolin 1; Synoviolin 1 isoform b; SYNOVIOLIN; SYVN1;;SYVN1
基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	Q86TM6.A synthesized peptide derived from human SYVN1

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:2000, ICC 1:50-1:200, FC 1:20-1:100, IP 1:20-1:50
蛋白分子量 (Molecular Weight)	Calculated MW: 68 kDa ; Observed MW: 65 kDa

研究背景 (Background)


产品名称: SYVN1/HRD1 Rabbit Monoclonal Antibody
产品货号: AMRe84525




E3 ubiquitin-protein ligase which accepts ubiquitin specifically from endoplasmic reticulum-associated UBC7 E2 ligase and transfers it to substrates, promoting their degradation.

研究领域 (Research Area)

图片 (Image Data)



All lanes use the Antibody at 1:2K dilution for 1 hour at room temperature.



Immunoprecipitate (IP) analysis using the Antibody at 1:50 dilution. (wb at 1:1K dilution)

注意事项 (Note)

For research use only .