

产品名称: Aminopeptidase A/CD249 Rabbit Monoclonal Antibody
产品货号: AMRe84534



产品概述 (Summary)

产品名称 (Production Name)	Aminopeptidase A/CD249 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide, 0.05% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	Aminopeptidase A/CD249
别名 (Alternative Names)	Aminopeptidase A; APA; bp1; CD249; EAP; Enpep; Gp160; Ly51;;Glutamyl aminopeptidase
基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	Q07075.A synthesized peptide derived from human Glutamyl aminopeptidase

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:2000,IHC 1:100-1:200
蛋白分子量 (Molecular Weight)	Calculated MW: 109 kDa ; Observed MW: 160 kDa

产品名称: Aminopeptidase A/CD249 Rabbit Monoclonal Antibody
产品货号: AMRe84534




研究背景 (Background)


Appears to have a role in the catabolic pathway of the renin-angiotensin system. Probably plays a role in regulating growth and differentiation of early B-lineage cells.

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of Aminopeptidase A / CD249 expression in human fetal lung lysate.



All lanes use the Antibody at 1:1K dilution for 1 hour at room temperature.

注意事项 (Note)

For research use only .