

产品名称: SRSF3 Rabbit Monoclonal Antibody
产品货号: AMRe84537



产品概述 (Summary)

产品名称 (Production Name)	SRSF3 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF,ICC,FC
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide,0.05% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	SRSF3
别名 (Alternative Names)	Splicing factor; SRp20; SRSF3; X16;;SRSF3
基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	P84103.A synthesized peptide derived from human SRSF3

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:2000,IHC 1:100-1:200,ICC/IF 1:50-1:200,ICC 1:50-1:200,FC 1:20-1:100
蛋白分子量 (Molecular Weight)	19 kDa

研究背景 (Background)


产品名称: SRSF3 Rabbit Monoclonal Antibody
产品货号: AMRe84537




Splicing factor, which binds the consensus motif 5'-C[ACU][AU]C[ACU][AC]C-3' within pre-mRNA and promotes specific exons inclusion during alternative splicing. Interaction with YTHDC1, a RNA-binding protein that recognizes and binds N6-methyladenosine (m6A)-containing RNAs, promotes recruitment of SRSF3 to its mRNA-binding elements adjacent to m6A sites within exons.

研究领域 (Research Area)


图片 (Image Data)



All lanes use the Antibody at 1:2K dilution for 1 hour at room temperature.



All lanes use the Antibody at 1:2K dilution for 1 hour at room temperature.



产品名称: SRSF3 Rabbit Monoclonal Antibody
产品货号: AMRe84537



Western blot analysis of SRSF3 expression in (1) Jurkat cell lysate; (2) Rat brain lysate.

注意事项 (Note)

For research use only .