

产品名称: GTPBP4 Rabbit Monoclonal Antibody
产品货号: AMRe84546



产品概述 (Summary)

产品名称 (Production Name)	GTPBP4 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,FC,IP
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide,0.05% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	GTPBP4
别名 (Alternative Names)	CRFG; GTP binding protein 4; Gtpbp4; NGB; NOG1;;GTPBP4
基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	Q9BZE4.A synthesized peptide derived from human GTPBP4

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:2000,IHC 1:100-1:200,FC 1:20-1:100,IP 1:20-1:50
蛋白分子量 (Molecular Weight)	Calculated MW: 74 kDa ; Observed MW: 80 kDa

研究背景 (Background)

Involved in the biogenesis of the 60S ribosomal subunit.Acts as a TP53 repressor, preventing TP53 stabilization and cell


产品名称: GTPBP4 Rabbit Monoclonal Antibody
产品货号: AMRe84546



cycle arrest.

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of GTPBP4 expression in Human fetal kidney lysate.

注意事项 (Note)

For research use only .