

产品名称: Calpain small subunit 1 Rabbit Monoclonal Antibody
产品货号: AMRe84549



产品概述 (Summary)

产品名称 (Production Name)	Calpain small subunit 1 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF,ICC,FC
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide,0.05% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	Calpain small subunit 1
别名 (Alternative Names)	Calpain4; CalpainS1; CANP; CANPS; CAPN4; CAPNS; Capns1; CDPS; CSS1;;Calpain small subunit 1
基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	P04632.A synthesized peptide derived from human Calpain small subunit 1

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:2000,IHC 1:100-1:200,ICC/IF 1:50-1:200,ICC 1:50-1:200,FC 1:20-1:100
蛋白分子量 (Molecular Weight)	28 kDa

产品名称: Calpain small subunit 1 Rabbit Monoclonal Antibody
产品货号: AMRe84549




研究背景 (Background)


Regulatory subunit of the calcium-regulated non-lysosomal thiol-protease which catalyzes limited proteolysis of substrates involved in cytoskeletal remodeling and signal transduction. Essential for embryonic development.

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of Calpain small subunit 1 expression in 293T cell lysate.



All lanes use the Antibody at 1:1K dilution for 1 hour at room temperature.

注意事项 (Note)

For research use only .