

产品名称: EWSR1 Rabbit Monoclonal Antibody
产品货号: AMRe84558



产品概述 (Summary)

产品名称 (Production Name)	EWSR1 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF,ICC
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide,0.05% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	EWSR1
别名 (Alternative Names)	EWS; Ewsr1; RNA binding protein EWS;;EWSR1
基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	Q01844.A synthesized peptide derived from human EWSR1

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:2000,IHC 1:100-1:200,ICC/IF 1:50-1:200,ICC 1:50-1:200
蛋白分子量 (Molecular Weight)	Calculated MW: 68 kDa ; Observed MW: 85 kDa

研究背景 (Background)

Might normally function as a transcriptional repressor. EWS-fusion-proteins (EFPS) may play a role in the tumorigenic


产品名称: EWSR1 Rabbit Monoclonal Antibody
产品货号: AMRe84558



process. They may disturb gene expression by mimicking, or interfering with the normal function of CTD-POLII within the transcription initiation complex. They may also contribute to an aberrant activation of the fusion protein target genes.

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of EWSR1 expression in K562 cell lysate.

注意事项 (Note)

For research use only .