

产品名称: Lactoferrin Rabbit Monoclonal Antibody
产品货号: AMRe84559



产品概述 (Summary)

产品名称 (Production Name)	Lactoferrin Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB, ICC, IP
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide, 0.05% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	Lactoferrin
别名 (Alternative Names)	GIG12; HLF2; Kaliocin 1; Lactoferrin; Lactotransferrin; LTF; Neutrophil lactoferrin; Talalactoferrin; Lactoferrin
基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	P02788.A synthesized peptide derived from human Lactoferrin

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:2000, ICC 1:50-1:200, IP 1:20-1:50
蛋白分子量 (Molecular Weight)	Calculated MW: 78 kDa ; Observed MW: 85 kDa

研究背景 (Background)


产品名称: Lactoferrin Rabbit Monoclonal Antibody
产品货号: AMRe84559



Transferrins are iron binding transport proteins which can bind two Fe^{3+} ions in association with the binding of an anion, usually bicarbonate.

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of Lactoferrin expression in HeLa cell lysate.

注意事项 (Note)

For research use only .