

产品名称: NRG3 Rabbit Monoclonal Antibody
产品货号: AMRe84561



产品概述 (Summary)

产品名称 (Production Name)	NRG3 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,FC
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide,0.05% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	NRG3
别名 (Alternative Names)	HRG3; Neuregulin 3; Nrg3; Pro neuregulin 3;;NRG3
基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	P56975.A synthesized peptide derived from human NRG3

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:2000,FC 1:20-1:100
蛋白分子量 (Molecular Weight)	Calculated MW: 78 kDa ; Observed MW: 56 kDa

研究背景 (Background)

Direct ligand for the ERBB4 tyrosine kinase receptor. Binding results in ligand-stimulated tyrosine phosphorylation and


产品名称: NRG3 Rabbit Monoclonal Antibody
产品货号: AMRe84561



activation of the receptor. Does not bind to the EGF receptor, ERBB2 or ERBB3 receptors. May be a survival factor for oligodendrocytes.

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of NRG3 expression in (1) SH-SY5Y cell lysate; (2) C6 cell lysate.

注意事项 (Note)

For research use only .