

产品名称: GSC Rabbit Monoclonal Antibody
产品货号: AMRe84565



产品概述 (Summary)

产品名称 (Production Name)	GSC Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IP
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide, 0.05% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	GSC
别名 (Alternative Names)	Goosecoid; GSC;;GSC
基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	P56915.A synthesized peptide derived from human GSC

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:2000, IP 1:20-1:50
蛋白分子量 (Molecular Weight)	28 kDa

研究背景 (Background)

Regulates chordin (CHRD). May play a role in spatial programming within discrete embryonic fields or lineage compartments


产品名称: GSC Rabbit Monoclonal Antibody
产品货号: AMRe84565




during organogenesis. In concert with NKX3-2, plays a role in defining the structural components of the middle ear; required for the development of the entire tympanic ring.

研究领域 (Research Area)


图片 (Image Data)



Western blot analysis of GSC expression in MCF7 cell lysate.



All lanes use the Antibody at 1:1K dilution for 1 hour at room temperature.



All lanes use the Antibody at 1:1k dilution for 1 hour at room temperature.

产品名称: GSC Rabbit Monoclonal Antibody
产品货号: AMRe84565



注意事项 (Note)

For research use only .