

产品名称: PSMA2 Rabbit Monoclonal Antibody
产品货号: AMRe84568



产品概述 (Summary)

产品名称 (Production Name)	PSMA2 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF,ICC,FC
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide,0.05% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	PSMA2
别名 (Alternative Names)	HC3; PMSA2; Proteasome component C3; PSC2; PSC3; psmA2;;PSMA2
基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	P25787.A synthesized peptide derived from human PSMA2

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:2000,IHC 1:100-1:200,ICC/IF 1:50-1:200,ICC 1:50-1:200,FC 1:20-1:100
蛋白分子量 (Molecular Weight)	Calculated MW: 26 kDa ; Observed MW: 25 kDa

研究背景 (Background)


产品名称: PSMA2 Rabbit Monoclonal Antibody
产品货号: AMRe84568



Component of the 20S core proteasome complex involved in the proteolytic degradation of most intracellular proteins. This complex plays numerous essential roles within the cell by associating with different regulatory particles. Associated with two 19S regulatory particles, forms the 26S proteasome and thus participates in the ATP-dependent degradation of ubiquitinated proteins.

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of PSMA2 expression in HepG2 cell lysate.

注意事项 (Note)

For research use only .