

产品名称: TIMM50 Rabbit Monoclonal Antibody
产品货号: AMRe84570



产品概述 (Summary)

产品名称 (Production Name)	TIMM50 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC
种属反应性 (Reactivity)	Human,Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide,0.05% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	TIMM50
别名 (Alternative Names)	TIM50; TIM50L; TIMM50;;TIMM50
基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	Q3ZCQ8.A synthesized peptide derived from human TIMM50

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:2000,IHC 1:100-1:200
蛋白分子量 (Molecular Weight)	Calculated MW: 40 kDa ; Observed MW: 36 kDa

研究背景 (Background)

Essential component of the TIM23 complex, a complex that mediates the translocation of transit peptide-containing


产品名称: TIMM50 Rabbit Monoclonal Antibody
产品货号: AMRe84570



proteins across the mitochondrial inner membrane. Has some phosphatase activity in vitro; however such activity may not be relevant in vivo.

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of TIMM50 expression in 293T cell lysate.

注意事项 (Note)

For research use only .