

产品名称: NEK7 Rabbit Monoclonal Antibody
产品货号: AMRe84571



产品概述 (Summary)

产品名称 (Production Name)	NEK7 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,FC
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide,0.05% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	NEK7
别名 (Alternative Names)	NEK7; NimA related protein kinase 7;;NEK7
基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	Q8TDX7.A synthesized peptide derived from human NEK7

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:2000,FC 1:20-1:100
蛋白分子量 (Molecular Weight)	35 kDa

研究背景 (Background)

Protein kinase which plays an important role in mitotic cell cycle progression.Required for microtubule nucleation activity


产品名称: NEK7 Rabbit Monoclonal Antibody
产品货号: AMRe84571




of the centrosome, robust mitotic spindle formation and cytokinesis.

研究领域 (Research Area)


图片 (Image Data)



All lanes use the Antibody at 1:2K dilution for 1 hour at room temperature.



All lanes use the Antibody at 1:2K dilution for 1 hour at room temperature.



Western blot analysis of NEK7 expression in Jurkat cell lysate.

产品名称: NEK7 Rabbit Monoclonal Antibody
产品货号: AMRe84571



注意事项 (Note)

For research use only .