

产品名称: RIP2 Rabbit Monoclonal Antibody
产品货号: AMRe84576



产品概述 (Summary)

产品名称 (Production Name)	RIP2 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide,0.05% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	RIP2 CARD carrying kinase; CARD3; CARDIAK; CCK; CLARP kinase; GIG30; Growth inhibiting gene 30; Receptor interacting protein 2; RICK; RIP 2; Ripk2; TNFRSF;;RIP 2
别名 (Alternative Names)	
基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	O43353.A synthesized peptide derived from human RIP 2

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:2000
蛋白分子量 (Molecular Weight)	61 kDa

产品名称: RIP2 Rabbit Monoclonal Antibody
产品货号: AMRe84576




研究背景 (Background)


Serine/threonine/tyrosine-protein kinase that plays an essential role in modulation of innate and adaptive immune responses. Acts as a key effector of NOD1 and NOD2 signaling pathways: upon activation by bacterial peptidoglycans, NOD1 and NOD2 oligomerize and recruit RIPK2 via CARD-CARD domains, leading to the formation of RIPK2 filaments.

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of RIP2 in K562 cell lysate.



All lanes use the Antibody at 1:2K dilution for 1 hour at room temperature.

注意事项 (Note)

For research use only .