

产品名称: ERAP1 Rabbit Monoclonal Antibody
产品货号: AMRe84579



产品概述 (Summary)

产品名称 (Production Name)	ERAP1 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide,0.05% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	ERAP1
别名 (Alternative Names)	ALAP; Aminopeptidase PILS; APPILS; Arts1; Endoplasmic reticulum aminopeptidase 1; ERAAP; ERAAP1; Erap1; PILSA; PILSAP;;ERAP1
基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	Q9NZ08.A synthesized peptide derived from human ERAP1

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:2000,IHC 1:100-1:200
蛋白分子量 (Molecular Weight)	Calculated MW: 107 kDa ; Observed MW: 100 kDa

研究背景 (Background)


产品名称: ERAP1 Rabbit Monoclonal Antibody
产品货号: AMRe84579




Aminopeptidase that plays a central role in peptide trimming, a step required for the generation of most HLA class I-binding peptides. Peptide trimming is essential to customize longer precursor peptides to fit them to the correct length required for presentation on MHC class I molecules. Strongly prefers substrates 9-16 residues long. Rapidly degrades 13-mer to a 9-mer and then stops. Preferentially hydrolyzes the residue Leu and peptides with a hydrophobic C-terminus, while it has weak activity toward peptides with charged C-terminus.

研究领域 (Research Area)

图片 (Image Data)




Western blot analysis of ERAP1 expression in K562 cell lysate.



All lanes use the Antibody at 1:2K dilution for 1 hour at room temperature.

产品名称: ERAP1 Rabbit Monoclonal Antibody
产品货号: AMRe84579



All lanes use the Antibody at 1:2K dilution for 1 hour at room temperature.

注意事项 (Note)

For research use only .