

产品名称: DBF4 Rabbit Monoclonal Antibody
产品货号: AMRe84580



产品概述 (Summary)

产品名称 (Production Name)	DBF4 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide, 0.05% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	DBF4
别名 (Alternative Names)	ASK; CHIF; DBF4; DBF4 zinc finger A; DBF4A; ZDBF1;;DBF4
基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	Q9UBU7.A synthesized peptide derived from human DBF4

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:2000
蛋白分子量 (Molecular Weight)	Calculated MW: 77 kDa ; Observed MW: 80 kDa

研究背景 (Background)

Regulatory subunit for CDC7 which activates its kinase activity thereby playing a central role in DNA replication and cell


产品名称: DBF4 Rabbit Monoclonal Antibody
产品货号: AMRe84580



proliferation. Required for progression of S phase. The complex CDC7-DBF4A selectively phosphorylates MCM2 subunit at 'Ser-40' and 'Ser-53' and then is involved in regulating the initiation of DNA replication during cell cycle.

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of DBF4 expression in Jurkat cell lysate.

注意事项 (Note)

For research use only .