

产品名称: DP1 Rabbit Monoclonal Antibody
产品货号: AMRe84581



产品概述 (Summary)

产品名称 (Production Name)	DP1 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,FC
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide,0.05% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	DP1
别名 (Alternative Names)	DP1; DRTF1; TFDP1;;DP1
基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	Q14186.A synthesized peptide derived from human DP1

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:2000,FC 1:20-1:100
蛋白分子量 (Molecular Weight)	Calculated MW: 45 kDa ; Observed MW: 50 kDa

研究背景 (Background)

Can stimulate E2F-dependent transcription. Binds DNA cooperatively with E2F family members through the E2 recognition


产品名称: DP1 Rabbit Monoclonal Antibody
产品货号: AMRe84581



site, 5'-TTTC[CG]CGC-3', found in the promoter region of a number of genes whose products are involved in cell cycle regulation or in DNA replication

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of DP1 expression in Raji cell lysate.

注意事项 (Note)

For research use only .