

产品名称: TMOD2 Rabbit Monoclonal Antibody
产品货号: AMRe84584



产品概述 (Summary)

产品名称 (Production Name)	TMOD2 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF,ICC,IP
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide,0.05% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	TMOD2
别名 (Alternative Names)	Tropomodulin 2; NTMOD; TMOD2;;Tropomodulin 2
基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	Q9NZR1.A synthesized peptide derived from human Tropomodulin 2

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:2000,IHC 1:100-1:200,ICC/IF 1:50-1:200,ICC 1:50-1:200,IP 1:20-1:50
蛋白分子量 (Molecular Weight)	40 kDa

研究背景 (Background)


产品名称: TMOD2 Rabbit Monoclonal Antibody
产品货号: AMRe84584



Blocks the elongation and depolymerization of the actin filaments at the pointed end. The Tmod/TM complex contributes to the formation of the short actin protofilament, which in turn defines the geometry of the membrane skeleton.

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of TMOD2 expression in SH-SY5Y cell lysate.

注意事项 (Note)

For research use only .