

产品名称: CDC123 Rabbit Monoclonal Antibody
产品货号: AMRe84585



产品概述 (Summary)

产品名称 (Production Name)	CDC123 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,FC,IP
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide,0.05% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	CDC123
别名 (Alternative Names)	C10orf7; CDC123; Cell division cycle 123; D123; HT 1080; PZ32;;CDC123
基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	O75794.A synthesized peptide derived from human CDC123

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:2000,FC 1:20-1:100,IP 1:20-1:50
蛋白分子量 (Molecular Weight)	39 kDa

研究背景 (Background)

ATP-dependent protein-folding chaperone for the eIF2 complex. Binds to the gamma subunit of the eIF2 complex which


产品名称: CDC123 Rabbit Monoclonal Antibody
产品货号: AMRe84585



allows the subunit to assemble with the alpha and beta subunits (By similarity).

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of CDC123 expression in HeLa cell lysate.

注意事项 (Note)

For research use only .