

产品名称: DOCK2 Rabbit Monoclonal Antibody
产品货号: AMRe84588



产品概述 (Summary)

产品名称 (Production Name)	DOCK2 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide, 0.05% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	DOCK2
别名 (Alternative Names)	DOCK 2; IMD40;;Dedicator of cytokinesis protein 2
基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	Q92608.A synthesized peptide derived from human Dedicator of cytokinesis protein 2

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:2000, IHC 1:100-1:200
蛋白分子量 (Molecular Weight)	212 kDa

研究背景 (Background)


产品名称: DOCK2 Rabbit Monoclonal Antibody
产品货号: AMRe84588



Involved in cytoskeletal rearrangements required for lymphocyte migration in response of chemokines. Activates RAC1 and RAC2, but not CDC42, by functioning as a guanine nucleotide exchange factor (GEF), which exchanges bound GDP for free GTP. May also participate in IL2 transcriptional activation via the activation of RAC2.

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of DOCK2 expression in Jurkat cell lysate.

注意事项 (Note)

For research use only .