

产品名称: STK39 Rabbit Monoclonal Antibody
产品货号: AMRe84589



产品概述 (Summary)

产品名称 (Production Name)	STK39 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,FC
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide, 0.05% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	STK39
别名 (Alternative Names)	DCHT; PASK; SPAK; Stk39;;STK39
基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	Q9UEW8.A synthesized peptide derived from human STK39

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:2000,IHC 1:100-1:200,FC 1:20-1:100
蛋白分子量 (Molecular Weight)	Calculated MW: 59 kDa ; Observed MW: 68 kDa

研究背景 (Background)

Effector serine/threonine-protein kinase component of the WNK-SPAK/OSR1 kinase cascade, which is involved in various


产品名称: STK39 Rabbit Monoclonal Antibody
产品货号: AMRe84589




processes, such as ion transport, response to hypertonic stress and blood pressure. Specifically recognizes and binds proteins with a RFXV motif.

研究领域 (Research Area)

图片 (Image Data)



All lanes use the Antibody at 1:2K dilution for 1 hour at room temperature.



All lanes use the Antibody at 1:2K dilution for 1 hour at room temperature.

注意事项 (Note)

For research use only .