

产品名称: CXCL5 Rabbit Monoclonal Antibody
产品货号: AMRe84591



产品概述 (Summary)

产品名称 (Production Name)	CXCL5 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB, ICC, IP
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide, 0.05% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	CXCL5
别名 (Alternative Names)	AMCFII; Cxcl5; ENA78; SCYB5;; CXCL5
基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	P42830.A synthesized peptide derived from human CXCL5

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:2000, ICC 1:50-1:200, IP 1:20-1:50
蛋白分子量 (Molecular Weight)	Calculated MW: 12 kDa ; Observed MW: 10 kDa

研究背景 (Background)

Involved in neutrophil activation. In vitro, ENA-78(8-78) and ENA-78(9-78) show a threefold higher chemotactic activity for


产品名称: CXCL5 Rabbit Monoclonal Antibody
产品货号: AMRe84591



neutrophil granulocytes.

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of CXCL5 expression in A549 cell lysate treated with TNF alpha and TPA .

注意事项 (Note)

For research use only .