

产品名称: EDG3 Rabbit Monoclonal Antibody  
产品货号: AMRe84592



## 产品概述 (Summary)

|                        |                                        |
|------------------------|----------------------------------------|
| 产品名称 (Production Name) | EDG3 Rabbit Monoclonal Antibody        |
| 描述 (Description)       | Recombinant rabbit monoclonal antibody |
| 宿主 (Host)              | Rabbit                                 |
| 应用 (Application)       | WB,IHC,ICC/IF,ICC                      |
| 种属反应性 (Reactivity)     | Human                                  |

## 产品性能 (Performance)

|                     |                                                                                             |
|---------------------|---------------------------------------------------------------------------------------------|
| 偶联物 (Conjugation)   | Unconjugated                                                                                |
| 修饰 (Modification)   | Unmodified                                                                                  |
| 同种型 (Isotype)       | IgG                                                                                         |
| 克隆 (Clonality)      | Monoclonal                                                                                  |
| 形式 (Form)           | Liquid                                                                                      |
| 存放说明 (Storage)      | Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.    |
| 储存溶液 (Buffer)       | Purified antibody in PBS with 0.05% sodium azide,0.05% protective protein and 50% glycerol. |
| 纯化方式 (Purification) | Affinity Purification                                                                       |

## 免疫原信息 (Immunogen)

|                        |                                                      |
|------------------------|------------------------------------------------------|
| 基因名 (Gene Name)        | EDG3                                                 |
| 别名 (Alternative Names) | EDG3; LPB3; S1P receptor 3; S1P3; S1pr3;;EDG3        |
| 基因 ID (Gene ID)        |                                                      |
| 蛋白 ID (SwissProt ID)   | Q99500.A synthesized peptide derived from human EDG3 |

## 产品应用 (Application)

|                          |                                                                   |
|--------------------------|-------------------------------------------------------------------|
| 稀释比 (Dilution Ratio)     | WB 1:1000-1:2000,IHC 1:100-1:200,ICC/IF 1:50-1:200,ICC 1:50-1:200 |
| 蛋白分子量 (Molecular Weight) | 42 kDa                                                            |

## 研究背景 (Background)

Receptor for the lysosphingolipid sphingosine 1-phosphate (S1P). S1P is a bioactive lysophospholipid that elicits diverse


**产品名称: EDG3 Rabbit Monoclonal Antibody**  
**产品货号: AMRe84592**



physiological effect on most types of cells and tissues. When expressed in rat HTC4 hepatoma cells, is capable of mediating S1P-induced cell proliferation and suppression of apoptosis.

## 研究领域 (Research Area)

## 图片 (Image Data)



Western blot analysis of EDG3 expression in HepG2 cell lysate.

## 注意事项 (Note)

For research use only .