

产品名称: HMGN2 Rabbit Monoclonal Antibody
产品货号: AMRe84595



产品概述 (Summary)

产品名称 (Production Name)	HMGN2 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide, 0.05% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	HMGN2
别名 (Alternative Names)	HMG17; HMGN2;;HMGN2
基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	P05204.A synthesized peptide derived from human HMGN2

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:2000, IHC 1:100-1:200
蛋白分子量 (Molecular Weight)	Calculated MW: 9 kDa ; Observed MW: 18 kDa

研究背景 (Background)

Binds to the inner side of the nucleosomal DNA thus altering the interaction between the DNA and the histone octamer.


产品名称: HMGN2 Rabbit Monoclonal Antibody
产品货号: AMRe84595



May be involved in the process which maintains transcribable genes in a unique chromatin conformation.

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of HMGN2 expression in K562 cell lysate.

注意事项 (Note)

For research use only .