

产品名称: IL11RA Rabbit Monoclonal Antibody
产品货号: AMRe84596



产品概述 (Summary)

产品名称 (Production Name)	IL11RA Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF,ICC,FC
种属反应性 (Reactivity)	Human,Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide,0.05% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	IL11RA
别名 (Alternative Names)	CRSDA; IL-11 receptor subunit alpha; IL-11R subunit alpha; IL11RA; Interleukin 11 receptor alpha; Interleukin 11 receptor alpha chain;;IL 11 receptor alpha
基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	Q14626.A synthesized peptide derived from human IL 11 receptor alpha

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:2000,IHC 1:100-1:200,ICC/IF 1:50-1:200,ICC 1:50-1:200,FC 1:20-1:100
蛋白分子量 (Molecular Weight)	45 kDa

产品名称: IL11RA Rabbit Monoclonal Antibody
产品货号: AMRe84596




研究背景 (Background)


Receptor for interleukin-11 (IL11). The receptor systems for IL6, LIF, OSM, CNTF, IL11 and CT1 can utilize IL6ST for initiating signal transmission. The IL11/IL11RA/IL6ST complex may be involved in the control of proliferation and/or differentiation of skeletogenic progenitor or other mesenchymal cells (Probable). Essential for the normal development of craniofacial bones and teeth. Restricts suture fusion and tooth number.

研究领域 (Research Area)

图片 (Image Data)




Western blot analysis of IL11RA expression in 293T cell lysate.



All lanes use the Antibody at 1:2K dilution for 1 hour at room temperature.

产品名称: IL11RA Rabbit Monoclonal Antibody
产品货号: AMRe84596



All lanes use the Antibody at 1:2K dilution for 1 hour at room temperature.

注意事项 (Note)

For research use only .