

产品名称: eIF3B Rabbit Monoclonal Antibody
产品货号: AMRe84598



产品概述 (Summary)

产品名称 (Production Name)	eIF3B Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB, ICC, FC, IP
种属反应性 (Reactivity)	Human, Mouse, Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide, 0.05% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	eIF3B
别名 (Alternative Names)	eIF3b; EIF3S9; hPrt1; PRT1;;eIF3B
基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	P55884.A synthesized peptide derived from human eIF3B

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:2000, ICC 1:50-1:200, FC 1:20-1:100, IP 1:20-1:50
蛋白分子量 (Molecular Weight)	Calculated MW: 92 kDa ; Observed MW: 116 kDa

研究背景 (Background)

Component of the eukaryotic translation initiation factor 3 (eIF-3) complex, which is required for several steps in the


产品名称: eIF3B Rabbit Monoclonal Antibody
产品货号: AMRe84598




initiation of protein synthesis.

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of eIF3B expression in (1) A431 cell lysate; (2) 293T cell lysate.



Immunoprecipitate (IP) analysis using the Antibody at 1:50 dilution. (wb at 1:3K dilution)

注意事项 (Note)

For research use only .