

产品名称: EEF1G Rabbit Monoclonal Antibody
产品货号: AMRe84599



产品概述 (Summary)

产品名称 (Production Name)	EEF1G Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF,ICC,FC,IP
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide,0.05% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	EEF1G
别名 (Alternative Names)	eEF 1B gamma; EEF1G; EF 1G; EF1 gamma; EF1G; GIG35;;eEF1B gamma
基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	P26641.A synthesized peptide derived from human eEF1B gamma

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:2000,IHC 1:100-1:200,ICC/IF 1:50-1:200,ICC 1:50-1:200,FC 1:20-1:100,IP 1:20-1:50
蛋白分子量 (Molecular Weight)	Calculated MW: 50 kDa ; Observed MW: 48 kDa

研究背景 (Background)


产品名称: EEF1G Rabbit Monoclonal Antibody
产品货号: AMRe84599




Probably plays a role in anchoring the complex to other cellular components.

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of EEF1G expression in (1) 293T cell lysate; (2) HepG2 cell lysate; (3) NIH/3T3 cell lysate.



Immunoprecipitate (IP) analysis using the Antibody at 1:50 dilution. (wb at 1:3K dilution)

注意事项 (Note)

For research use only .