

产品名称: RPL7A Rabbit Monoclonal Antibody
产品货号: AMRe84607



产品概述 (Summary)

产品名称 (Production Name)	RPL7A Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide,0.05% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	RPL7A
别名 (Alternative Names)	PLA X polypeptide; RP23-414L19.4; RPL7A; SURF3; surfet 3; TRUP;;RPL7A
基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	P62424.A synthesized peptide derived from human RPL7A

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:2000,IHC 1:100-1:200
蛋白分子量 (Molecular Weight)	30 kDa

研究背景 (Background)

Component of the large ribosomal subunit. The ribosome is a large ribonucleoprotein complex responsible for the


产品名称: RPL7A Rabbit Monoclonal Antibody
产品货号: AMRe84607



synthesis of proteins in the cell.

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of RPL7A expression in (1) HeLa cell lysate; (2) Raw264.7 cell lysate.

注意事项 (Note)

For research use only .